



Emerging Technologies in Sample Analysis

4 April 2002

New England Bioterrorism Preparedness Workshop
MIT Lincoln Laboratory

Bernadette Johnson
(781)981-1902
bernadette@ll.mit.edu

MS-15453


This work was sponsored under Air Force contract F19628-00-C-0002. The views expressed are those of the Author and do not reflect official policy or position of the United States Government.

50
years



Report Documentation Page

Report Date 04APR2002	Report Type N/A	Dates Covered (from... to) 03APR2002 - 04APR2002
Title and Subtitle Emerging Technologies in Sample Analysis	Contract Number F19628-00-C-0002	
	Grant Number	
	Program Element Number	
Author(s) Johnson, Bernadette	Project Number	
	Task Number	
	Work Unit Number	
Performing Organization Name(s) and Address(es) MIT Lincoln Laboratory	Performing Organization Report Number	
Sponsoring/Monitoring Agency Name(s) and Address(es) Air Force ESC/XPK (Richard Axtell) Hanscom AFB, MA 01731	Sponsor/Monitor's Acronym(s)	
	Sponsor/Monitor's Report Number(s)	
Distribution/Availability Statement Approved for public release, distribution unlimited		
Supplementary Notes Workshop paper from the New England Bioterrorism Preparedness Workshop held 3-4 april 2002 at MIT Lincoln Laboratory, Lexington, MA, The original document contains color images.		
Abstract		
Subject Terms		
Report Classification unclassified	Classification of this page unclassified	
Classification of Abstract unclassified	Limitation of Abstract SAR	
Number of Pages 19		



Outline

- **Current techniques in sample analysis**
 - Clinical (subject of yesterday's talk)
 - Environmental
- **Challenges associated with environmental sampling**
- **Examples of technologies in use and in development**




CDC's Sample Analysis Guidelines


(example: *B. Anthracis*)

- **Persons suspected of exposure/infection**
 - Cultures of blood and spinal fluid
 - Cultures of tissues or fluids from affected areas
 - Microscopic examination
 - PCR
 - Nasal swab (occasionally for exposure, but not for diagnosis)
 - Antibody testing (exposure, not validated for diagnosis)

- **Environmental contamination**
 - Cultures of air samples, surface swabs, suspicious powders
 - Microscopic examination of suspect material
 - Evaluation of growth properties of suspect agent
 - PCR
 - DFA (direct fluorescent assay) to detect key bacterial proteins
 - Specialized tests, such as immunoassays (SMART)



How Do These Techniques Compare?





Examples of In-use and **Developmental** Immunoassay Devices




← 4.5 in →



**Ticket cartridges and reader
for lateral-flow immunoassay
in Joint Biological Point Detection System
(JBPDS)**




**Response Equipment Co.
Bio-HAZ Biodetector**



Dendrimer-Based Alert Ticket (ARL)


Disposable Test Strip



Prototype Design




**Upconverting
Phosphors (SRI)**




Features of Immunoassay Analysis

- Can be used on environmental samples with little or no preparation
- Readout is fast (~ 15 minutes) and simple (colorimetric or fluorimetric)
- Sensitivity modest (~10,000 - 100,000 particles)
 - Depends on antibody-antigen binding affinity and readout scheme
- Specificity reasonably good
 - Depends on antibody construct and antigen specificity
- Current IAs are not multiplexed; development of protein microarrays may lead to sensitive, multi-assay analysis tools





Examples of Existing Protein Microarrays



Phylos (2000 element)




Ciphergen (multiple classes of proteins)




Zyomyx (10,000/cm²)

- Protein microarray technology development driven by drug screening and disease-marker investigations
 - Diagnostics (clinical and environmental) still developmental



Developmental Antibody-Based Sensor: CANARY


Concept



Prototype microcentrifuge device



Tests Against Killed Tularemia (Collab. with NMRC)




Status of B-Cell Lines

Complete

FMDV
VEE
Vibrio cholera
Orthopox viruses
Yersinia pestis
Brucella spp
Francisella tularensis

In development

Coxiella burnetti
Bacillus anthracis
E. coli O157:H7



PCR-Based Analysis Tools

- **Systems being developed (and deployed) that provide agent ID within 30 minutes of introduction of prepared sample**

Semi-automated field-portable
PCR devices



RAPID - Idaho Technologies




***SmartCycler
XC System - Cepheid***

Example of handheld PCR
device




***HANAA - Handheld Nucleic
Acid Analyzer, developed by
LLNL, Cepheid, and ETG, Inc.***

- **Challenge remains in automating sample preparation and analysis**
 - Pathogen cells or spores must be ruptured to liberate the DNA/RNA
 - DNA/RNA must be separated from protein debris/environmental impurities




Overview of Sample Preparation






Target Concentration: Affinity Magnet Protocol





DNA Purification: Simple Nucleic Acid Prep (SNAP)

- Chemically treated paper is the key component of SNAP
- Lyses cells, binds PCR-assay inhibitors, and purifies DNA
- **Advantages:**
 - Fast and easy (1/5th the time of other published protocols)
 - Water is only added reagent (no phenol, chloroform, or alcohol)
 - Lightweight, compact, enables archiving
 - On-site fixation: preserves DNA & kills pathogenic organisms






Lincoln Interim Nucleic-acid Kit (LINK)




(Developed in response to October 2001 events)

- **LINK as a solution:**
 - Incorporates SNAP paper but in a more user-friendly format
 - Faster processing than basic SNAP
 - Easier to sample, handle, and process
 - Enables on-site fixation
 - Outside can be decontaminated
 - 6 minute processing time
 - Single-step processing
 - Results equal to or better than basic SNAP






How to Use LINK



1) Apply sample
Sit for 5 minutes


2) Process in **one step**


3) Remove DNA
Total time ~6 minutes!



LINK Cartridge Works with Varied Samples


- LINK detection from:
 - **Portal Shield** air-to-liquid samples seeded with vegetative bacteria
 - Untreated domestic sewage (Boston) seeded with vegetative bacteria
 - Paper, envelopes, skin seeded with bacterial spores
 - Air impaction with dry bacterial spores





What About DNA Microarrays?

- DNA Microarray: Any 2D or 3D substrate having many ($\sim 10^2$ - 10^5) different nucleic-acid capture sites (probes)
- Can identify both strain and drug resistance of pathogens
- Can offer highly multiplexed assay capability






Pathogen Identification via DNA Microarray

- Detect small amounts (<100 copies per ml) of pathogen-specific nucleic acids in environmental sample
- Arrays might provide log orders more information than current PCR-based approaches (e.g. TaqMan)
- Challenges for diagnostic applications:
 - Never demonstrated for environmental (or clinical) samples
 - Amplification may be necessary before micro-array assay
 - Sample preparation required (as in PCR techniques)



Assay Times for Current and Emerging PCR/DNA Systems





Summary

- **Environmental sample analysis parallels methodology developed for clinical sampling**
 - Immunoassays for rapid estimate of exposure (not yet CDC authorized)
 - PCR techniques being deployed in some laboratories to provide strain specificity and drug resistance
 - Culture still used to provide “gold standard” for pathogen ID
- **New technology developments could greatly increase the speed, sensitivity, and multiplicity of environmental assays**
 - Protein microarrays could offer highly multiplexed, rapid ID capability on collected samples
 - DNA microarrays could offer hundreds to thousands of pathogen tests on single-chip format, provided sample preparation can be made compatible